

## Aerobic Culture

## Anaerobic Culture

Aerobic Culture				Anaerobic Culture			
Order	Source/Site	Container	Transport Temp	Order	Source/Site	Container	Transport Temp
Culture, Body Fluid <b>LABBFCULT</b>	Any fluid Identify source/site in EPIC	Sterile syringe with rubber cap <b>OR</b> Sterile urine cup	Refrigerated	Anaerobic Culture <b>LAB6080</b>	Any Identify source/site in EPIC	T588 (ACTI glass tube with black cap)	Ambient and <b>upright (Must be placed on top of the media inside the glass tube. Do not push it into the media)</b>
Culture, Wound <b>LAB6007</b>	Any swab Identify source/site in EPIC	E-swab with white cap ( <b>preferred</b> ) <b>OR</b> Red double swab culturette	Refrigerated				
Culture, Tissue <b>LAB6903</b>	Tissue or Bone Identify source/site in EPIC	Sterile urine cup with enough saline to cover tissue (0.5mL)	Refrigerated				
Please see RMC catalog for other culture types or for additional information							