





Lab Reference Grid – CT/NG & TV media

EPIC Test Code	Swab/Container	Image	Source	Performing Lab
LAB6059 (CTGCM)	COBAS Dual Swab Kit		Woven Swab= vaginal collections Flocked Swab= endocervical collections, please use woven swab as cleaning swab	TTMC
LAB6059 (CTGCM)	COBAS Urine Kit		First catch “dirty” urine 10-50mL Note: must wait 1 hour if patient recently voided or collected a clean catch sample	TTMC
LABTVMDU	Sterile UA container		First catch “dirty” urine 20-50mL Note: must wait 1 hour if patient recently voided or collected a clean catch sample	Waconia
LABCGRNA			Rectal or Pharyngeal	Mayo