


Lab Reference Grid – CT/NG Media

EPIC Test Code	Swab	Image	Source	Performing Lab
LAB6059 (CTGCM)	COBAS Dual Swab Kit		Woven swab = vaginal collections Flocked swab = endocervical collections, please use woven swab as cleaning swab	TTMC
LAB6059 (CTGCM)	COBAS Urine Kit		Urine	TTMC
LABCGRNA	Aptima Collection Multitest Swab or Aptima Unisex Swab Kit		Rectal or Pharyngeal	Mayo
LABMTRNA (Male only)	Aptima Urine		Useful for: Trichomonas vaginalis in urine and male patient specimens Urine	Mayo